Semi-Micro Osmometer K-7400

Total Osmolality Determination of Biological Fluids

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The new **Semi-Micro Osmometer K-7400** measures the **freezing point depression** of liquids and it is based on the successful predecessor model Type Dig. L

Principle

The principle of total osmolality is to determine the osmolality of body fluids such as blood (serum) and urine. This is an essential routine process in clinics, institutes, and universities. Osmolality describes the concentration of osmotic effective particles in solutions, independent from type, composition or electrical charge. The osmolality refers to the mass of the solution, that means to 1 kg pure water.

Advantages

- Simple handling
- 2 to 4 point calibration
- Short measurement time (2-3 min.)
- Reasonable price
- Fully automated measurement
- Automatic storing of measurement values
- Air cooling
- Cooling process is thermistor controlled
- Cooling chamber will not freeze

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Medical Applications

- Poisoning
- Monitoring kidney functions
- Diabetes Insipidus
- Shock (trauma)
- Clearance control
- Lactation control
- Dialyse control
- Control of infusion solutions in pharmacies
- Control of iso, hyper, and hypotonic solutions
- Blood control and control of other transfusion supplies
- Control of osmotic pressure in cell culture media for biotechnology and genetic technology

Operation

150 µl of an aqueous solution or a biological fluid are super cooled below the melting point. Then a stirring needle (vibrator) automatically initiates crystallization. The crystallization warmth causes a temperature rises up to the melting point (below 0°C). The temperature difference to 0°C measured is proportional to the number of free particles in the solution. The measurement value in mOsmol is shown in the display.

Technical Data

Measurement Accuracy:	error ≤ 1%
Measurement Precision :	RSD < 1%
Calibration :	2 calibration points necessary up to 4 calibration points accepted
Sample Volume :	150 µl
Measuring Time :	2 minutes
Guaranteed Range:	0 – 2000 mOsmol
Possible Range:	0 – 3500 mOsmol
RS-232 Interface Connection to Printer A0591 or to PC	
Power Supply :	90-260 V, 47 - 63Hz, 70 VA
Dimensions (W/H/D):	160 x 182 x 340 mm
Printer:	76.8 x 77.4 x 39.3 mm
Weight :	5.0 kg
Printer :	0.15 kg

Ordering Information

Order No.

A3707	Semi-Micro Osmometer K-7400
A3707-1	Semi-Micro Osmometer K-7400 with printer
A3709	Semi-Micro Osmometer K-7400 with Measuring Head for plastic vials
A3709-1	Semi-Micro Osmometer K-7400 with printer and with Measuring Head for plastic vials

Accessories

A0591	Printer
A0840	Measuring Head for glass vials
A0840-1	Measuring Head for plastic vials
Y1241	Calibration Solution 400mOsmol
Y1240	Calibration Solution 300mOsmol
A0913	Measuring Vials (glass, 150 µl), 1000 pcs
A0720	Measuring Vials (plastic, 150 µl), 1000 pcs

Technical data are subject to change without notice.

Visit www.knauer.net for details on complete HPLC systems, HPLC columns, and osmometers.

KNAUER · ASI · Advanced Scientific Instruments
Wissenschaftliche Gerätebau
Dr. Ing. Herbert Knauer GmbH
Hegauer Weg 38, 14163 Berlin, Germany
Telephone: +49 (0)30 / 80 97 27-0
Telefax: +49 (0)30 / 8 01 50 10
E-Mail: info@knauer.net
Internet: www.knauer.net

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